

Pelvi Center Contra-Indications

Device:	QRS 1010 PelviCenter V 1.2
Release company:	QRS International AG
Date of release:	1.3.2013

1. Contra-indications Absolutely

Pregnancy

If a patient is pregnant or the patient plans to become pregnant during the treatment period, it is forbidden to participate the treatment.

Metal Implants between knee and neck

- Female hormone spiral,
- hip or knee replacements,
- screws, nails,
- etc.

Patients often mean that only materials which can become magnetic will be a no go. It is important to understand that **each kind of metal is a no go**.

Electronically controlled Implants between knee and neck

- Insulin pump
- Pace maker
- Neuron stimulator
- etc.

Strong magnetic fields have a potential to disturb and interference electronic devices positioned in the field.. Using devices of the above list or similar **is a not go**

Non Metal Implants between knee and neck

- meshes
- etc.

The PelviCenter will have a very strong effect on your muscles in the whole urogenital track. The arranged organs are affected in a manner not predictable. And so maybe implants like meshes will do.

It is strongly recommended that patients will do a discussion with their responsible doctors, to clarify, whether a treatment with the PelviCenter is possible or not.

The operator is not able to decide.

If the patient is affected by the contra-indication and has no written confirmation from his doctor, PelviCenter treatment is not permitted.

Surgical interventions

If patients would have had surgery interventions in body areas affected by the magnetic field, it is strongly recommended that patients will do a discussion with their responsible doctors, to clarify, whether a treatment with the PelviCenter is possible or not.

The operator is not able to estimate the potential risk the patient would have.

If the patient is affected by the contra-indication and has no written confirmation from his doctor, PelviCenter treatment is not permitted.

Severe cardiac arrhythmias

If patients would have had severe cardiac arrhythmias, it is strongly recommended that patients will do a discussion with their responsible doctors, to clarify, whether a treatment with the PelviCenter is possible or not.

The operator is not able to estimate the potential risk the patient would have.

If the patient is affected by the contra-indication and has no written confirmation from his doctor, PelviCenter treatment is not permitted.

Epilepsy

If patients would have epilepsy, it is strongly recommended that patients will do a discussion with their responsible doctors, to clarify, whether a treatment with the PelviCenter is possible or not.

The operator is not able to estimate the potential risk the patient would have.

If the patient is affected by the contra-indication and has no written confirmation from his doctor, PelviCenter treatment is not permitted.

2. Contra-indications Relatively

The following contra-indications are no absolutely no go's.

But we know from our Patients that patients under these circumstances might have a lot of uncomfortable feelings during the treatment.

- Menstrual
- Acute urinary tract infection
- Painful hemorrhoids
- Febrile infection